

Swing-out rotor 4758



illustrated with
carriers 1369-92











Swing-out rotor, 2-place















illustrated with carriers 4745
and lids 4627 with
bio-containment ⁵⁾ and
removal frame 4626

∠ 90° n = 5,100 min⁻¹ max. RCF 3,926

Cat. No. (without carriers) 4728

capacity in ml	1.5	2.0	1	3	4
Ø x L in mm	11x38		6x45	10x60	12x60
Cat. No.	2078	0536	Rhesus	-	-
					
carrier Cat. No. 1366					
Cat. No.	5277	1357	1327	1326	
boring Ø x L in mm	11.5x38	6.5x23	10.5x23	12.5x44	
tubes per rotor	72	240	96		
max. RCF ²⁾	2,470	2,454	2,438		
radius in mm	153	152	151		
run-up in sec	19				
run-down in sec, braked	16				
temperature in °C ¹⁾	0				

T x B x H in mm	86x128x15	86x128x17.5	86x128x22	86x128x44.5	86x128x46	86x128x83	59x84x11	82x124x20	-
capacity in ml	-								0.2
Cat. No.	MTP		CP	DWP	MS	QP	Microtest plate (Terasaki)	PCR plate, 96 wells	PCR strips
lid with bio-containment ⁵⁾ Cat. No. 4627									
rotor 4728 carrier Cat. No. 4745									
Cat. No.	4626							4626 + 1485	
boring Ø x L in mm	-								
plates / strips per rotor	12	10	8	2		4	2	24x8	
max. RCF ²⁾	3,926								
radius in mm	135								
run-up in sec	67								
run-down in sec, braked	41								
temperature in °C ¹⁾	+4								

¹⁾ Packed in units of 10 pieces.

²⁾ Please load only the middle row.