






































Swing-out rotor, 8-place	Swing-out rotor, 8-place
	
\sphericalangle 90° n = 4,000 min ⁻¹ max. RCF 2,415	\sphericalangle 45° n = 4,000 min ⁻¹ max. RCF 2,469
Cat. No. (without carriers) 1611	Cat. No. (with carriers) 1617

capacity in ml	5	6	7	10	15	2.6–3.4	2.7–3	4–5.5	4.5–5	7.5–8.5	1.6–5	4–7	4–7	8.5–10	
Ø x L in mm	12x75	12x82	12x100	17x70	13x100	17x100	13x65	11x66	15x75	11x92	15x92	13x75	13x100	16x75	16x100
Cat. No.	0553 ²⁾	0501 ²⁾	0578 ²⁾	2079	-	0518 ²⁾	blood collection / urine tubes								
 rotor Cat. No. 1611															
															
Cat. No.	1131	1643	1132	1643	1644	1131	1132	1643	1644	1131	1643	1132	1643	1132	1644
boring Ø x L in mm	13x53	13.2x81	17.5x53	13.2x81	17.5x81	13x53	17.5x53	13.2x81	17.5x81	13x53	13.2x81	17.5x53	13.2x81	17.5x53	17.5x81
tubes per rotor	8														
max. RCF ²⁾	1,914	2,415	1,914	2,415	2,415	1,914	1,914	2,415	2,415	1,914	2,415	1,914	2,415	1,914	2,415
radius in mm	107	135	107	135	135	107	107	135	135	107	135	107	135	107	135
run-up in sec	22														
run-down in sec, braked	25														

capacity in ml	15	50
Ø x L in mm	17x120	29x115
Cat. No.	0509	0513
 rotor Cat. No. 1617		
		
Cat. No.	1462-A	-
boring Ø x L in mm	17x84	30x94.5
tubes per rotor	8	
max. RCF ²⁾	2,469	
radius in mm	138	
run-up in sec	22	
run-down in sec, braked	25	