

Swing-out rotor, 8-place

Swing-out rotor, 12-place



∠ 45°

n = 5,000 min⁻¹
max. RCF 3,857





Cat. No. (with carriers) **1617**
























∠ 55°
(with carriers 1127)
∠ 60°
(with carriers 1122)
∠ 80°
(with carriers 1621)

n = 5,000 min⁻¹
max. RCF 4,193

Cat. No. (without carriers) **1628**

capacity in ml	15	50
∅ x L in mm	17x120	29x115
Cat. No.	0509	0513
 rotor Cat. No. 1617		
		
Cat. No.	1462-A	-
boring ∅ x L in mm	17x84	30x94.5
tubes per rotor	8	
max. RCF ²⁾	3,857	
radius in mm	138	
run-up in sec	20	
run-down in sec, braked	19	
temperature in °C ¹⁾	-10	

capacity in ml	5	10	15	2.6–3.4	2.7–3	4–5.5	7.5–8.5	1.6–5	4–7	8.5–10
∅ x L in mm	12x75	17x70	17x100	13x65	11x66	15x75	15x92	13x75	16x75	16x100
Cat. No.	0553²⁾	2079	0518²⁾	blood collection / urine tubes						
 rotor Cat. No. 1628										
										
Cat. No.	1127	1122	1621	1127	1122	1621	1127	1122	1621	
boring ∅ x L in mm	13.2x53	17.5x53	17.5x79	13.2x53	17.5x53	17.5x79	13.2x53	17.5x53	17.5x79	
tubes per rotor	12									
max. RCF ²⁾	3,494	3,522	4,193	3,494	3,522	4,193	3,494	3,522	4,193	
radius in mm	125	126	150	125	126	150	125	126	150	
run-up in sec	16									
run-down in sec, braked	16									
temperature in °C ¹⁾	-15	-10	-15	-15	-10	-15	-15	-10	-15	