

Angle rotor, 14-place

Angle rotor, 6-place, for Schlenk tubes



∠ 36°

n = 9,000/10,000 min⁻¹
max. RCF 14,127/17,441

with bio-containment⁵⁾

Cat. No. 5655



∠ 45°

n = 2,000 min⁻¹
max. RCF 805

Cat. No. 5616³⁴⁾

capacity in ml	15	25	7.5–8.5	9–10	10	8.5–10	15	50	30
∅ x L in mm	17x100	24x100	15x92	16x92	15x102	16x100	17x120	29x115	26x95
Cat. No.	0518²⁾	0519²⁾	blood/urine tubes				0509	0513	0545
lid with bio-containment ⁵⁾ incl.									
rotor Cat. No. 5655									
Cat. No.	1473	1472	1473				1474	-	1472
boring ∅ x L in mm	17.5x91.5	26x85	17.5x91.5				17x107	30x100	26x85
tubes per rotor	14								
max. RCF ²⁾ classic/cooled	13,493/ 16,658	13,131/ 16,211	13,493/ 16,658				13,765/ 16,994	14,127/ 17,441	13,131/ 16,211
radius in mm	149	145	149				152	156	145
run-up in sec classic/cooled	124/ 145								
run-down in sec, braked classic/cooled	94/ 108								
temperature in °C ¹⁾	+16								

capacity in ml	25	50
∅ x L in mm	24x146.5	38x148.5
Cat. No.	0532	0533
 rotor Cat. No. 5616³⁴⁾		
Cat. No.	4317	-
boring ∅ x L in mm	25x86	38.5x89.9
tubes per rotor	6	
max. RCF	783	805
radius in mm	175	180
run-up in sec	17	
run-down in sec, braked	21	
temperature in °C ¹⁾	-14	

¹³⁾ The screw cap for tube No. 0538 is available with Cat. No. 0539. Dimensions with cap 38x110 mm.

³⁴⁾ In the heatable models 5660-50, 5660-51, 5660-70, 5670-50 and 5675-50, rotor 5616 may only be used up to a maximum of +40 °C.