


































| Angle rotor, 6-place | Angle rotor, 14-place |
|---|--|
|  |  |
| \sphericalangle 25° | \sphericalangle 36° |
| n = 9,500 min ⁻¹ max. RCF 14,025 | n = 10,000 min ⁻¹ max. RCF 17,441 |
| Cat. No. 4266 | Cat. No. 4248 |

| | | | | | | | | | |
|--|---|---|---|---|---|--|---|---|---|
| capacity in ml | 15 | 25 | 94 | 10 | 30 | 50 | 85 | 94 | 250 |
| Ø x L in mm | 17x100 | 24x100 | 38x102 | 16x80 | 26x95 | 29x107 | 38x106 | 38x102 | 61x122 |
| Cat. No. | 0518²⁾ | 0519²⁾ | 0549 | - | 0545 | 0546 | 0547 | 0538¹³⁾ | 5127²⁸⁾ |
| lid with bio-containment ⁵⁾ incl. |  |  |  |  |  |  |  |  |  |
|  rotor Cat. No. 4266 |  |  |  |  |  |  |  |  | |
| Cat. No. | 5646 | 5642 | 5644 | 5641 | 5642 | 5643 | 5644 | | - |
| boring Ø x L in mm | 17.6x83 | 26x80 | 38.6x88 | 16.6x70 | 26x80 | 29x90 | 38.6x88 | | 61.5x109 |
| tubes per rotor | 42 | 18 | 6 | 48 | 18 | 6 | | | |
| max. RCF ²⁾ | 13,319 | 12,915 | 12,310 | 13,420 | 12,915 | 12,108 | 12,310 | | 14,025 |
| radius in mm | 132 | 128 | 122 | 133 | 128 | 120 | 122 | | 139 |
| run-up in sec | 82 | | | | | | | | |
| run-down in sec, braked | 96 | | | | | | | | |
| temperature in °C ¹⁾ | +2 | | | | | | | | |

| | | | | | | | | | |
|--|---|---|---|---|---|---|---|---|---|
| capacity in ml | 15 | 25 | 7.5-8.5 | 9-10 | 10 | 8.5-10 | 15 | 50 | 30 |
| Ø x L in mm | 17x100 | 24x100 | 15x92 | 16x92 | 15x102 | 16x100 | 17x120 | 29x115 | 26x95 |
| Cat. No. | 0518²⁾ | 0519²⁾ | blood collection/urine tubes | | | | 0509 | 0513 | 0545 |
| lid with bio-containment ⁵⁾ incl. |  |  |  |  |  |  |  |  |  |
|  rotor Cat. No. 4248 |  |  |  | | | |  | |  |
| Cat. No. | 1473 | 1472 | 1473 | | | | 1474 | - | 1472 |
| boring Ø x L in mm | 17.5x91.5 | 26x85 | 17.5x91.5 | | | | 17x107 | 30x100 | 26x85 |
| tubes per rotor | 14 | | | | | | | | |
| max. RCF ²⁾ | 16,658 | 16,211 | 16,658 | | | | 16,994 | 17,441 | 16,211 |
| radius in mm | 149 | 145 | 149 | | | | 152 | 156 | 145 |
| run-up in sec | 103 | | | | | | | | |
| run-down in sec, braked | 106 | | | | | | | | |
| temperature in °C ¹⁾ | +6 | | | | | | | | |

¹³⁾ The screw cap for tube No. 0538 is available with Cat. No. 0539. Dimensions with cap 38x110 mm.

²⁸⁾ At temperatures of over +40 °C and/or when not filled to capacity, bottles may warp during centrifugation.